

Cataract Suture Needle Ranges

Needle	Suture Material	Suture Length	Suture Sizes		
			10-0	9-0	8-0
Advanced Micro Point Spatula AS - 140 - 6 Diameter - 0.15mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	15cm (6")	6401N	6491N	
		30cm (12")	6402N	6492N	
Advanced Micro Point Spatula AS - 160 - 6 Diameter - 0.15mm Length - 5.6mm ½ Circle - 160°	aurolon Nylon Black Monofilament	15cm (6")	6601N	6691N	
		30cm (12")	6602N	6692N	
Advanced Micro Point Spatula AS - 140 - 8 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	40cm (16")	8401N	8491N	8481N
		40cm (16")	8402N	8492N	8482N
Advanced Micro Point Spatula AS - 140 - 8 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	aurosilk Silk Black Braided Silicone Coated	15cm (6")			8481S
		30cm (12")			8482S
Advanced Micro Point Spatula AS - 140 - 8 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	aurosilk Silk Blue Virgin	40cm (16")			8482VSB*
Advanced Micro Point Spatula AS - 140 - 0606 Diameter - 0.15mm Length - 6.0mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Monofilament	10cm (4")	A0101*		

Scleral Fixation Surgery Suture Needle

Needle	Suture Material	Suture Length	Suture Sizes		
			10-0	9-0	7-0
Advanced Micro Point Spatula AS - 000 - 6 Diameter - 0.15mm Length - 16.0mm Straight Needle	aurolene Polypropylene Blue Monofilament	20cm (8")	6002PP		
Advanced Micro Point Spatula AS - 000 - 6 Diameter - 0.15mm Length - 16.0mm Straight Needle	aurolene Polypropylene Blue Monofilament	20cm (8")		6092PP	
Advanced Micro Point Spatula AS - 000 - 6, AS - 160 - 6 Diameter - 0.15mm Length - 16mm, 5.6mm Straight Needle, ½ Circle - 160°	aurolene Polypropylene Blue Monofilament	30cm (12")	60602PP		
Advanced Micro Point Spatula AS - 000 - 6, AS - 160 - 6 Diameter - 0.15mm Length - 5.6mm, 16.0mm Straight Needle, ½ Circle - 160°	aurolene Polypropylene Blue Monofilament	30cm (12")		60692PP	
Advanced Micro Point Spatula AS - 090 - 0616 Diameter - 0.15mm Length - 16mm, Straight Needle	aurolene Polypropylene Blue Monofilament	30cm (12")			6072PP
Advanced Micro Point Spatula AS - 000 - 6, AS-160-6 Diameter - 0.15mm Length - 16mm, 16.2mm, Straight Needle, 1/2 Circle 160°	aurolene Polypropylene Blue Monofilament	30cm (12")			60674PP

Iris Fixation Surgery Suture Needle

Advanced Micro Point Spatula AS - 090 - 6 Diameter - 0.15mm Length - 16 mm 1/4 Circle - 90°	aurolene Polypropylene Blue Monofilament	20cm (8")	6902PP*				
--	---	-----------	---------	--	--	--	--

Corneal Suture Needle Ranges

Needle	Suture Material	Suture Length	Suture Sizes					
			10-0	9-0	8-0	7-0	6-0	4-0
Advanced Micro Point Spatula AS - 140 - 6 Diameter - 0.15mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	15cm (6")	6401N	6491N				
		30cm (12")	6402N	6492N				
Advanced Micro Point Spatula AS - 160 - 6 Diameter - 0.15mm Length - 5.6mm ½ Circle - 1	aurolon Nylon Black Monofilament	15cm (6")	6601N	6691N				
		30cm (12")	6602N	6692N				
Advanced Micro Point Spatula AS - 140 - 8 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	40cm (16")	8401N	8491N				
		40cm (16")	8402N	8492N				
Advanced Micro Point Spatula AS - 135 - 1108 Diameter - 0.28mm Length - 8.0mm 3/8 Circle - 135°	aurosilk Silk Black Braided Silicone Coated	76cm (30")					B6105	
		76cm (30")					B6106	
Conventional Cutting CC - 135 - 1716 Diameter - 0.43mm Length - 15.9mm 3/8 Circle - 135°	aurosilk Silk Black Braided Silicone Coated	76cm (30")					B4401	
Advanced Micro Point Spatula AS - 090-1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	PolyCryl Poly Glycolic Acid Violet Braided	45cm (18")					A6108*	
Advanced Micro Point Spatula AS - 140 - 0806 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30cm (12")			A8104*			
Advanced Micro Point Spatula AS - 140 - 0606 Diameter - 0.15mm Length - 6.0mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Monofilament	10cm (4")	A0101*					
Advanced Micro Point Spatula AS - 140 - 0907 Diameter - 0.230mm Length - 6.5mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30cm (12")					A7102*	

Keratoplasty Surgery Suture Needle

Advanced Micro Point Spatula AS - 180 - 6 Diameter - 0.15mm Length - 5.6mm ½ Circle - 180°	aurolon Nylon Black Monofilament	30cm (12")	6802N*				
--	---	------------	--------	--	--	--	--






Orbit & Oculoplast Suture Needle Ranges

Needle	Suture Material	Suture Length	Suture Sizes				
			8-0	7-0	6-0	5-0	4-0
Advanced Micro Point Spatula AS - 135 - 1108 Diameter - 0.28mm Length - 8.0mm 3/8 Circle - 135°	aurosilk Silk Black Braided Silicone Coated	76cm (30")			B6105		
		76cm (30")			B6106		
Conventional Cutting CC - 135 - 1716 Diameter - 0.43mm Length - 15.9mm 3/8 Circle - 135°	aurosilk Silk Black Braided Silicone Coated	76cm (30")					B4401
Reverse Cutting RC - 130 - 1513 Diameter - 0.38mm Length - 13.1mm 3/8 Circle - 130°	aurolene Polypropylene Blue Monofilament	70cm (28")					P5301
Advanced Micro Point Spatula AS - 140 - 0806 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30cm (12")	A8104*				
Advanced Micro Point Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	PolyCryl Poly Glycolic Acid Violet Braided	45cm (18")				A6108*	
Advanced Micro Point Spatula AS - 140 - 0907 Diameter - 0.230mm Length - 6.5mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30 cm (12")		A7102*			
Taper Point TP - 180 - 1716 Diameter - 0.430mm Length - 16mm ½ Circle - 180°	aurosilk Silk Black Braided Silicone Coated	90 cm (36")					B4203*
Advanced Micro Point Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	PolyCryl Poly Glycolic Acid Violet Braided	45cm (18")				A5108*	












Squint & Membranectomy Suture Needle Ranges

Needle	Suture Material	Suture Length	Suture Sizes	
			6-0	5-0
Advanced Side Cutting Spatula AS - 90 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	PolyCryl Poly Glycolic Acid Violet Braided	45cm (18")	A6108*	
Reverse Cutting RC - 130 - 1513 Diameter - 0.38mm Length - 13.1mm 3/8 Circle - 130°	aurolene Polypropylene Blue Monofilament	70cm (28")		P5301

Retinal Suture Needle Ranges

Needle	Suture Material	Suture Length	Suture Sizes			
			7-0	6-0	5-0	4-0
Advanced Micropoint Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	PolyCryl Poly Glycolic Acid Violet Braided	45cm (18") 		A6108 *		
Micro Point Spatula AS - 090 - 1308 Diameter - 0.33mm Length - 8.0mm 1/4 Circle - 90°	aurobond Polyester Green Braided	76cm (30") 			E5101	
		76cm (30") 			E5102	
Micro Point Spatula AS - 090 - 1708 Diameter - 0.43mm Length - 7.9mm 1/4 Circle - 90°	aurobond Polyester Green Braided	76cm (30") 				E4101
Advanced Micropoint Spatula AS - 140 - 0907 Diameter - 0.230mm Length - 6.5mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30 cm (12") 	A7102 *			

Glaucoma Surgery Suture Needle

Needle	Suture Material	Suture Length	Suture Sizes			
			10-0	8-0	7-0	3-0
Advanced Micro Point Spatula AS - 140 - 6 Diameter - 0.15mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	15cm (6") 	6401N			
		30cm (12") 	6402N			
Advanced Micro Point Spatula AS - 160 - 6 Diameter - 0.15mm Length - 5.6mm 1/2 Circle - 160°	aurolon Nylon Black Monofilament	15cm (6") 	6601N			
		30cm (12") 	6602N			
Advanced Micro Point Spatula AS - 180 - 6 Diameter - 0.15mm Length - 5.6mm 1/2 Circle - 180°	aurolon Nylon Black Monofilament	30cm (12") 	6802N *			
Advanced Micro Point Spatula AS - 140 - 8 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	aurolon Nylon Black Monofilament	40cm (16") 	8401N	8481N		
		40cm (16") 	8402N	8482N		
Taper Point TP - 130 - 0607 Diameter - 0.15mm Length - 6.7mm 3/8 Circle - 130°	aurolon Nylon Black Monofilament	30cm (12") 	N0204			
Advanced Micro Point Spatula AS - 140 - 0806 Diameter - 0.20mm Length - 6.0mm 3/8 Circle - 140°	PolyCryl Poly Glycolic Acid Violet Braided	30cm (12") 		A8104 *		
Conventional Cutting CC - 135 - 1716 Diameter - 0.43mm Length - 15.9mm, 3/8 Circle - 135°	Supralon Nylon Black Multifilament	10cm (4") 				N3403
Reverse Cutting RC - 140 - 0908 Diameter - 0.23mm Length - 7.8mm, 3/8 Circle - 140°	aurosilk Silk Black Braided Silicone Coated	45cm (18") 			B7302	

Product Features

- Perfectly biocompatible products
- Advanced micro point spatula/taper/conventional and reverse cutting edges for minimal penetration force
- Burr free needle surface for lesser tissue trauma
- Special suture needle attachment process for maximum strength
- High suture tensile strength for continuous workability and knot security
- Optimum needle cutting edges for drag free passage through the toughest tissue
- Unique needle body geometry for excellent grip stability
- High needle bending resistance for retaining original curvature
- Homogenous suture surface for less tissue distortion

Flexible product characteristics to meet the surgeon's need in product ranges

Quality Policy

We shall strive for excellence in all our activities. We shall manufacture superior quality products and deliver on time with world class service at competitive prices to achieve total customer satisfaction, through continuous improvement by team work.

"Nylon, Silk & 2-0 to 8-0 PGA products are USFDA approved"




AUROLAB

Information published in this catalogue is subject to change without notification
 No.1, Sivagangai Main Road, Veerapanjan, Madurai - 625020,
 Tamilnadu, India. ☎ +91 94892 12354 | +91 452 3596100
 ✉ info@aurolab.com 🌐 www.aurolab.com



Issue No : 20 - 02/24

Ophthalmic Suture Needles



AUROLAB
SEE NEW POSSIBILITIES